

Alternative Western Blotting Protocol

(Produces less background than 4% BSA blocking)

- 1) Transfer gel to immobilization P membrane.
- 2) Block for 1-3 hrs in:
 - 2 ml 10X Western buffer
 - 16 ml dH₂O
 - 2 ml Horse SerumFilter sterilize or add 0.2% NaN₃ and save.
- 3) Incubate with primary antibody (1:500 dilution) overnight at 4°C in:
 - 4 ml 10X Western buffer
 - 28 ml dH₂O
 - 8 ml Horse SerumFilter sterilize or add 0.2% NaN₃ and save.
- 4) Wash (2 rinses, 1x15', 2x10') in:
 - 20 ml 10X Western buffer
 - 175 ml dH₂O
 - 5 ml 20% Tween 20
- 5) Incubate with secondary antibody (1:10,000) for 40-60min at r.t. in:
 - 4 ml 10X Western buffer
 - 28 ml dH₂O
 - 8 ml Horse Serum
- 6) Wash (2 rinses, 1x15', 2x10') in:
 - 24 ml 10X Western buffer
 - 84 ml dH₂O
 - 12 ml 5% Tween 20
- 7) Develop blot using Amersham ECL reagents as described by manufacturer.

10X Western buffer (1 L)

1.5 M NaCl (87.66g)

150 mM Tris pH 7.4 (150 ml of 1 M or 18.2g)